

## **APPLICATION NOTE MW-104**

## **CLOSED VESSEL DIGESTION**

Sample Tea Leaf

**Type of Vessel** HP/VHP Vessels and TFM Liners

Number of Vessels Up to 10 Reagents 8 mL HNO<sub>3</sub>

Sample Weight<sup>1</sup>: 0.2 grams Tea leaf

## **Procedure**

• Place 0.2 grams of the sample into the liner of each vessel

- Add small amount of water to wet the sample
- Add 8 ml of concentrated HNO<sub>3</sub>, to each vessel inside the fume hood and allow the sample to outgas before sealing the vessels.
- Pre-heat on the thermal plate at about 50 °C for 10 min.
- Load the vessel carousel into the microwave digestion oven and run digestion program as follows:

## **Digestion Program**

Step	Time/s	Starting Temp/C°	Ending Temp/C°
1	300	37	100
2	300	100	180
3	600	180	185
4	0	185	0

- After the digestion, the solution should be clear in slight green color.
- Add 1.5 ml of H<sub>2</sub>O<sub>2</sub> to each vessel and shake it well until it become colorless

NOTE: This procedure is a reference point for sample digestion using the Aurora Microwave digestion system and may need to be modified or changed to obtain the required results on your subject.

www.aurora-instr.com

