

Rb⁺ Efflux Assay for Assessment of Non-Selective Cation Channel Activities

Kun Liu, Manoj Samuel, Richard K. Harrison, and Jeff W. Paslay
Screening Sciences, Wyeth Research, Collegeville, Pennsylvania.

ABSTRACT

Transient receptor potential (TRP) channels have been found to play important roles in cellular physiology and hold promise as therapeutic targets. These channels activate in response to a variety of chemical or physical stimuli and conduct non-selective cation currents (NSCC). Due to their unique activation properties, application of automated electrophysiology to measure the channel activity has been difficult. Using HEK293 cells stably expressing human TRP channels, hTRPC6 and hTRPA1, we developed and validated a high-throughput Rb⁺ efflux assay for NSCC channels. The assay was performed in cell-based 96-well format. A significant increase in Rb⁺ efflux can be detected upon channel activation by specific agonists, confirming that both TRPC6 and TRPA1 channels are permeable to Rb⁺ ions. The agonists induced Rb⁺ efflux can be blocked by known channel blockers and selected compounds from our high-throughput screening (HTS) hits. The assay is suitable for HTS with Z' factors of 0.53 and above. We also tested the Ca²⁺ effect on channel activities in this assay. Both TRPC6 and TRPA1 channels were found to be inhibited by increasing the concentration of Ca²⁺ in the assay buffer. However, Ca²⁺ significantly reduced the potency of allyl isothiocyanate (AITC) on TRPA1 but did not affect the potency of carbachol on TRPC6. Using this assay for secondary confirmation screen, we successfully identified and confirmed the positive hits as TRPC6 inhibitors.

INTRODUCTION

Ion channels represent a large family of the promising therapeutic targets. During the last decade, an increasing number of ion channel proteins have been identified to play important roles in physiological and/or pathophysiological processes.

Among them, transient receptor potential (TRP) channels are the most promising candidates for drug intervention of disease processes and thus are actively pursued in pharmaceutical discovery.¹⁻³ TRP channels are not conventional voltage- or ligand-gated channels. They are gated by a variety of intracellular or extracellular signals. Upon activation, TRP channels carry non-selective cation currents (NSCC).^{4,5} One of the most interesting TRP channels is TRPC6, which belongs to canonical TRP channel subfamily. TRPC6 is physiologically activated through the activation of G-protein-coupled receptors linked with G_q pathway and can be directly opened by diacyl glycerol (DAG), a phospholipase C (PLC) hydrolyzed product.⁶ Recently, TRPC6 has been found to play an important role in cellular calcium homeostasis, especially during cell differentiation and proliferation. A more impressive finding is the up-regulation of this channel in some of the diseases related to tissue overgrowth, suggesting that it might be a potential target for diseases such as cardiac hypertrophy and cancers.^{3,7-13}

Despite the significant improvement made in ion channel assay technologies, it still lacks robust high-throughput assay for TRP channels, especially for TRPC6 channel. Currently available high-throughput assay is mostly fluorescence-based using Fluorometric Image Plate Reader platform (FLIPR). Although this assay is applicable for some of the TRP channels, the fluorescence labeling-related issue inevitably generates higher nonspecific hit rate resulting from the interaction with autofluorescent or quenching compounds. To address this issue, it is necessary to have a non-fluorescence-based assay for either primary or follow-up confirmation screen. The automated electrophysiology platforms such as IonWorks and Q-Patch have quickly evolved in recent years and have been successfully applied to screening of some ion channels including both voltage- and ligand-gated classes.^{14,15} However, due to the unique activation properties of TRP channels, this technology has not been successfully applied for these channels. Searching for a high-throughput, robust, and cost-efficient assay for screening on TRP channels is absolutely

ABBREVIATIONS: AITC, allyl isothiocyanate; CCH, carbachol; DAG, diacyl glycerol; FLIPR, Fluorometric Image Plate Reader; HTS, high-throughput screening; NMDG, N-methyl-D-glucamine; NSCC, non-selective cation currents; TRP, transient receptor potential.

necessary for identification of novel compounds against these potential drug targets.

Rubidium is a metal element with the physicochemical properties very similar to potassium. Rb^+ is often used as a tracer ion to assess potassium channel activities. The radioactive $^{86}\text{Rb}^+$ was broadly used for measuring K^+ channel function and/or K^+ transportation across biological membranes.^{16,17} In recent years, a nonradioactive label-free technology for Rb^+ measurement, flame atomic absorption spectrometry has been developed in automated format and successfully applied to the high-throughput assays for K^+ channel targets.^{18–20} The assays have been demonstrated sensitive, robust, and cost efficient. However, it is still unknown whether the assay can be utilized for NSCC channels. Using stable cell lines expressing TRP channels, we have demonstrated in the present study that Rb^+ efflux assay can be a robust high-throughput screening (HTS)-suitable assay format for functional assessment of NSCC channel activities. The assay has been successfully applied to our efforts in searching TRP channel modulators.

MATERIALS AND METHODS

Media and Chemicals

Cell culture media. For HEK293-hTRPC6, Dulbecco's modified Eagle's medium (DMEM) with 10% fetal bovine serum (FBS), 1% penicillin/streptomycin, and 150 $\mu\text{g}/\text{mL}$ hygromycin (Invitrogen, Carlsbad, CA) was used in cell culture. For HEK293-hTRPA1, DMEM/glutamax high glucose with 10% dialyzed FBS, 1% penicillin/streptomycin, 1% sodium pyruvate, 1% MEM NEAA, and 400 $\mu\text{g}/\text{mL}$ hygromycin was used.

Patch clamp pipette solution. The patch clamp pipette solution contains (in mM) 140 CsCl, 4 MgCl_2 , 10 EGTA, and 10 HEPES (pH 7.2).

Patch clamp external solution. The patch clamp external solution contains (in mM) 135 NaCl, 2 KCl, 1 CaCl_2 , 1.17 MgSO_4 , and 10 HEPES (pH 7.4). For testing Rb^+ permeability of the channel, NaCl was replaced by RbCl. For testing cation permeability, NaCl was replaced by *N*-methyl-D-glucamine (NMDG).

Rb^+ efflux assay buffer. For Rb^+ efflux assay, the assay buffer contains (in mM) 140 NaCl, 0.15 CaCl_2 , 3.3 KH_2PO_4 , 0.8 K_2HPO_4 , 1.2 MgCl_2 , 10 D-glucose, and 20 HEPES (pH 7.4).

Rb^+ efflux lysis buffer. 1% Triton® X-100 solution was used as lysis buffer. HC030031 was purchased from TOCRIS Bioscience (Ellisville, MO). *N*-(*p*-amylcinnamoyl)anthranilic acid (ACA) was purchased from EMD Chemicals Inc. (Gibbstown, NJ). All other chemicals and compounds were purchased from Sigma-Aldrich (St. Louis, MO). Compounds were prepared in dimethyl sulfoxide (DMSO) for stock and freshly diluted into the buffer with DMSO concentrations below 0.2%.

Stable Cell Line Generation

TRPC6 stable cell line was generated using HEK293 cells transfected with the human TRPC6 channel expression construct. hTRPC6 full-length DNA (accession number: BC093660; Thermo Fisher Scientific, Huntsville, AL) was subcloned into pCDNA-hygro vector (Invitrogen, Carlsbad, CA). After sequencing confirmation, the construct DNA was amplified and used to transfect HEK293 cells. HEK293 cells were seeded into a 10-cm culture dish at a density of $3.5\text{--}4 \times 10^6$ cells. After overnight incubation, cells were transfected with DNA construct using Lipofectamine™ 2000 (Invitrogen, Carlsbad, CA) according to the manufacturer's protocol. Five hours later, transfection medium was replaced and the cells were incubated overnight at 37°C in 95/5% air/ CO_2 . The transfected cells were then harvested and replated into a 96-well cell culture plate at a density of 2–5 cells/well in the medium containing 150 $\mu\text{g}/\text{mL}$ hygromycin. After 48 h of incubation, colonies were observed. The cells were then individually picked and scaled up for functional determination of channel expression using FLIPR membrane potential assay. TRPA1 stable cells were obtained from our colleagues in Department of Neuroscience, Wyeth Research.

Patch Clamp Recording

Cells were dissociated using trypsin (0.05 mg/mL) and harvested into 15-mL centrifuge tube. The freshly isolated cells were used within 4 h after preparation. High-resolution whole-cell current recordings were acquired with a computer-based patch clamp amplifier system (EPC-10; HEKA, Lambrecht, Germany). Patch pipettes had resistance between 2 and 5 megohm after filling with the standard pipette solution. After establishing whole-cell configuration, the membrane potential was held at -20 mV. A 1-s ramp voltage protocol from -100 to 100 mV for TRPC6 or from -60 to 60 mV for TRPA1 was applied repetitively to evoke channel currents in an interval of 5 s. Liquid junction potential, capacitive currents, and series resistance were determined and corrected using the automatic compensation of EPC-10. Compounds were delivered through a gravity-driven perfusion system. The current amplitudes at ± 60 mV were used for determining compound effect.

Rb^+ Efflux Assay

For Rb^+ efflux assay, cells were plated in 96-well poly-D-lysine-coated cell plates (BD, Franklin Lakes, NJ) at 50,000/well. After 48–72 h of incubation at 37°C in 95/5% air/ CO_2 , the culture medium was removed. One hundred microliters of Rb^+ loading buffer containing 5 mM RbCl was added into the wells. For constructing Rb^+ loading curve, cellular Rb^+ content was determined immediately after addition of Rb^+ loading buffer in every 8 wells

subsequently at each 20 min. For rest of the experiments to determine Rb⁺ efflux, 3 h of incubation in Rb⁺ loading buffer was applied for Rb⁺ loading process. After Rb⁺ loading, the loading buffer was removed and cells were washed twice with 200 μ L of plain assay buffer and then subjected to the compound treatment. The buffer containing various concentrations of the testing compounds was added 10 min before addition of the agonists at the indicated concentrations, followed by 30 min of agonist incubation for TRPC6 and 10 min for TRPA1. To measure Rb⁺ content, the supernatant was transferred into an empty plate and the cells were lysed in 150 μ L of the lysis buffer. Ten minutes later, 100 μ L of samples from both supernatant and cell lysis plates were used to determine the Rb⁺ content using ICR8000 automated Rb⁺ detection spectrometry (Aurora Biomed, Vancouver, Canada). Total Rb⁺ was calculated from addition of the Rb⁺ content in supernatant and cell lysis. To present the data as percent cellular Rb⁺ content, the following equation was applied:

$$\text{Percent cellular Rb}^+ \text{ content} = \frac{\text{Rb}^+ \text{ content in cell lysis}}{\text{total Rb}^+ \text{ content}} \times 100$$

Starting from Rb⁺ loading process, all experiments were performed at room temperature (22°C \pm 1°C).

Data Analysis

The Rb⁺ loading curves were fitted to single exponential growth equation as follows:

$$C_{\text{Rb}^+} = C_T \times (1 - \text{EXP}[-t/\tau])$$

where C_{Rb^+} is the Rb⁺ content at the time point t ; C_T is the maximal Rb⁺ content that can be loaded; τ is the time constant for the loading. For Rb⁺ efflux curves, the following equation of single exponential decay was used to get fitting parameters:

$$C_{\text{Rb}^+} = C_T \times \text{EXP}(-t/\tau)$$

The dose-dependent response was fitted using normal Hill equation to generate IC₅₀ (reduced signal) or EC₅₀ (increased signal). All data analyses were performed using the software Origin® 6.1 (Origin Lab Corporation, Northampton, MA). Data are presented in mean \pm SD.

RESULTS AND DISCUSSION

HEK293-hTRPC6 stable cell line and HEK293-hTRPA1 stable cell line are used in this study. These 2 channels are members of TRP family and share some common features such as nonselective permeability but have the least amino acid identity within the family (<20%). To confirm the channel expression and validate the cell lines for assay development, we performed patch

clamp recording to characterize the channel currents. As shown in Figure 1A, HEK293-hTRPC6 cell had a small current (<200 pA at \pm 60 mV) at resting state. The application of carbachol (CCH), an M cholinergic receptor agonist, dramatically increased both inward and outward currents. The CCH effect appeared quickly




Fig. 1. Patch clamp recording of TRPC6 stable cell line. **(A)** The current amplitudes evoked by voltage at -60 and 60 mV are plotted against time. After the formation of whole-cell configuration, a ramp voltage from -100 to 100 mV for 1 s was applied repetitively in an interval of 5 s. Carbachol (CCH) was added into the recording chamber at the given time indicated with the solid bar. **(B)** The current amplitudes evoked by the ramp voltage are plotted as a function of the membrane voltage. The representative current traces were recorded in the absence (control) and presence of CCH (CCH 100 μ M), in the buffer where Na⁺ was replaced by Rb⁺ (Rb⁺ + CCH 100 μ M), and in the buffer where Na⁺ was replaced by *N*-methyl-D-glucamine (NMDG + CCH 100 μ M).

but declined immediately after reaching the peak value. The declining process was slow and incomplete during 8-min treatment despite only small current retained. The reapplication of CCH evoked a small current increase, likely due to the depletion of the receptor channel coupling components. *Figure 1B* shows the voltage dependence of the currents obtained by the ramp voltage protocol. The CCH-induced current reversed at -3.65 ± 1.2 mV ($n = 8$), indicating that it is a nonselective current. The current amplitude was greater in the outward direction than in the inward direction and a slight dip was seen at the voltage from 0 to 30 mV, resembling the characteristic hallmarks of TRPC6 channels. Although the current was small without CCH, it had indeed the same voltage-dependent kinetics as the CCH-induced current, suggesting that the channel had basal opening. To examine Rb⁺ permeability of the channel, the external Na⁺ was replaced by the equal concentration of Rb⁺. No noticeable shift of reversal potential was observed, indicating that the channel was equally permeable to Na⁺ and Rb⁺. However, when Na⁺ was replaced by NMDG, the reversal potential was shifted to -62.1 ± 5.4 mV ($n = 5$), indicating that the channel was only permeable to cations. We were able to record the CCH-induced current with the amplitude ranging from 0.3 to 2 nA at 60 mV from 67% of the cells (8 of 12).

For HEK293-TRPA1 cell line, we did the same experiment to characterize the current induced by known TRPA1 agonist, allyl isothiocyanate (AITC). As shown in *Figure 2A*, HEK293-TRPA1 cell had little current before the addition of AITC. AITC increased both inward and outward currents evoked by ± 60 mV membrane voltage. This process was relatively slow but the current was sustained. No spontaneous decline was observed. The increase in extracellular Ca²⁺ further potentiated the current and then suddenly decreased and abolished it. This phenomenon has been previously described as Ca²⁺ potentiation and inactivation.²¹⁻²³ *Figure 2B* shows the voltage dependence of the current derived by the ramp voltage protocol. The AITC-induced current reversed at near 0 mV, indicating that this channel conducts a nonselective current. The channel was highly potentiated by Ca²⁺ before being completely inactivated, resembling the unique property of TRPA1 channel. Similar to TRPC6 channel, TRPA1 was also equally permeable to Na⁺ and Rb⁺ with no noticeable shift in reversal potential when Na⁺ was replaced by Rb⁺. NMDG significantly shifted the reversal potential to the left side, indicating the cationic permeability of the channel. The current evoked by 20 μ M AITC appeared in every cell tested with the amplitude in the range of 0.4 to 3 nA at 60 mV ($n = 12$).

Next, we started the Rb⁺ efflux assay development by examining the Rb⁺ loading process. In the presence of 5 mM of Rb⁺ in the loading buffer, the increase in intracellular Rb⁺ in both cell




Fig. 2. Patch clamp recording of TRPA1 stable cell line. **(A)** The current amplitudes evoked by voltage at -60 and 60 mV are plotted against time. The ramp voltage for TRPA1 is from -60 to 60 mV for 1 s. **(B)** The current amplitudes evoked by the ramp voltage are plotted as a function of the membrane voltage. The representative current traces were recorded in the absence (control) and presence of allyl isothiocyanate (AITC; AITC 20 μ M), in the presence of AITC and Ca²⁺ (AITC 20 μ M + Ca²⁺ 2 mM), in the buffer where Na⁺ was replaced by Rb⁺ (Rb⁺ + AITC 20 μ M), and in the buffer where Na⁺ was replaced by *N*-methyl-D-glucamine (NMDG + AITC 20 μ M).

lines followed the single exponential growth equation (*Fig. 3A* and *3B*). The time constants for HEK293-hTRPC6 and HEK293-TRPA1 are 45.5 and 64.9 min, respectively. The faster loading in TRPC6 cells may reflect the basal opening of the TRPC6 channels as described earlier.

Based on the loading kinetics, we chose the 3-h loading period for the Rb⁺ efflux determination. We then examined the basal Rb⁺ efflux in 1 h after Rb⁺ loading and the effect of channel agonists, CCH and AITC for TRPC6 and TRPA1, respectively, on Rb⁺ efflux. As shown in *Figure 4A*, under control condition, cellular




Fig. 3. Rb⁺ loading time course. Cellular Rb⁺ content was measured every 20 min after the cells were incubated with 5 mM RbCl. For each graph, the data were collected from 8 wells for each data point in a 96-well plate. The solid line represents the fitting curve using single exponential growth equation. **(A)** TRPC6 cell loading curve. **(B)** TRPA1 cell loading curve.

Rb⁺ content decreased time dependently with a time constant of 13.8 min. CCH induced a significant increase in the rate of reduction. Fitting to the single exponential decay curve resulted in a time constant of 11.7 min. Fifty minutes after treatment of CCH, the cellular Rb⁺ content decreased to 45% of the control. *Figure 4B* shows the Rb⁺ efflux in TRPA1 cells in the absence and presence of AITC. Similar to TRPC6 cells, the cellular Rb⁺ content in




Fig. 4. Time-dependent Rb⁺ efflux induced by the transient receptor potential (TRP) channel agonists. Immediately after washout of the loading buffer, Rb⁺ content was determined at each 5 min in 1 h except the first interval of 10 min. After the first 2 time point determinations, the agonists were added into the wells for agonist-treated group and the results were plotted together with the control group without agonist addition. For each data point, $n = 8$. The time course for both control and agonist-treated groups were fitted with single exponential decay equation and plotted in smooth solid lines. **(A)** Rb⁺ efflux in TRPC6 cells in the absence (open circle) and presence (solid circle) of carbachol (CCH). **(B)** Rb⁺ efflux in TRPA1 cells in the absence (open circle) and presence (solid circle) of AITC.

TRPA1 declined following the single exponential decay, but the rate was slower with the time constant of 25.1 min, consistent with the observation that TRPA1 showed much less basal current than TRPC6. AITC at 50 μM induced a quick decline with the time constant of 3.9 min. After 30 min treatment with AITC, the cellular Rb^+ content decreased to 8% of the control. The faster and more complete Rb^+ efflux in TRPA1 cells was in agreement with the results acquired in patch clamp recording.

The dose-dependent effects of the agonists were examined in both cell lines. To test the specificity of the response, we cross




Fig. 5. Dose-dependent effects of the TRP channel agonists. The cellular Rb^+ content was plotted against the concentration of both channel agonists. The data are fitted with Hill equation. (A) TRPC6 cells. (B) TRPA1 cells.

examined the specific agonists in both cell lines. As shown in *Figure 5A*, CCH dose dependently decreased the cellular Rb^+ content in TRPC6 cells with the EC_{50} at 3.6 μM . AITC had no effect in the same concentration range. In contrast, AITC evoked a dose-dependent reduction of the cellular Rb^+ content in TRPA1 cells with the EC_{50} at 6.2 μM but had no effect on TRPC6 cells (*Fig. 5B*).




Fig. 6. Effect of Ca^{2+} on Rb^+ efflux induced by channel agonists. The experiment was performed in the presence of the different concentrations of Ca^{2+} as indicated and the dose-response curves were established as described in *Figure 5*, except that the data presented here are percent cellular Rb^+ contents (see Materials and Methods). (A) TRPC6 cells. (B) TRPA1 cells.




Fig. 7. Effects of the TRPC6 channel blockers on CCH-induced Rb⁺ efflux in TRPC6 cells. **(A)** Dose-dependent inhibition of 100 μM CCH-induced Rb⁺ efflux by known TRPC6 channel blockers, SKF-96365 and *N*-(*p*-amylcinnamoyl)anthranilic acid (ACA). **(B)** Dose-dependent inhibition by 2 proprietary compounds. **(C)** A plate view of the raw data from a 96-well plate treated with control and 100 μM ACA. Shown on right are the bar-graphed statistic results for control wells ($n = 42$) and the blocker wells ($n = 42$).

These results further demonstrate that the responses induced by the agonists are specific to the corresponding channel activities.

We next tested the Ca²⁺ effect on the Rb⁺ efflux since both TRPC6 and TRPA1 have been reported to be regulated by calcium ions. We increased the Ca²⁺ concentration to 1.5 and 5 mM from 0.15 mM that we used in the earlier experiments. As shown in *Figure 6A*, increasing Ca²⁺ in assay buffer for TRPC6 shifted the dose-response curve up without changing the potency of CCH effect, suggesting that Ca²⁺ inhibits the channel activity regardless of CCH effect. But, increasing Ca²⁺ reduced the potency of AITC in TRPA1 cells (*Fig. 6B*), indicating that Ca²⁺ may affect AITC-binding properties.

Using 100 μM of CCH as agonist in the presence of 0.15 mM of Ca²⁺, we examined the effect of the known TRPC6 blockers, SKF-96365, and ACA on TRPC6 cells. As shown in *Figure 7A*, both blockers dose dependently inhibited the response induced by CCH. The IC₅₀s for these 2 compounds are 19.5 and 8.3 μM, respectively, comparable with the IC₅₀s acquired in FLIPR membrane potential assay (6.2 and 14 μM, respectively, unpublished data). The calculated Z' factors using maximal inhibition values are 0.67 and 0.55 for the 2 compounds, indicating that the assay is suitable for HTS (>0.5). *Figure 7B* shows the dose-response curves for our lead compounds derived from the primary HTS screen. These 2 compounds displayed potent dose-dependent inhibitions with IC₅₀s at 0.34 and 0.72 μM, nearly 20-folds greater than the potencies of the known blockers. These values are also in accordance with the IC₅₀s from FLIPR assay (0.52 and 0.72 μM, respectively, unpublished data). *Figure 7C* shows the raw data from a representative 96-well plate. The plate was treated with the buffer control and 100 μM ACA before addition of 100 μM CCH. The percent cellular Rb⁺ contents are 22.1 ± 3.2 and 55.4 ± 2.1, respectively. The Z' factor is 0.53, indicating that the assay is robust for HTS.

For TRPA1 assay, we examined the effect of a relatively selective channel blocker, HC030031, on AITC-induced Rb⁺ efflux.²⁴ Because Ca²⁺ could significantly reduce AITC effect, we performed this assay in the presence of 0.15 and 1.5 mM of Ca²⁺ for comparison. Using 50 μM of AITC as agonist, HC030031 caused a dose-dependent increase in cellular Rb⁺ content under both Ca²⁺ conditions (*Fig. 8A*). However, increase in Ca²⁺ concentration remarkably shifted the dose-response curve to the left side. The IC₅₀ was changed from 60 to 4.1 μM. Meanwhile, the assay window was significantly reduced from 51% to 27% in cellular Rb⁺ content due to a significant reduction in AITC response in the presence of 1.5 mM of Ca²⁺. Then we




Fig. 8. Effect of HC030031 on Rb⁺ efflux induced by AITC in TRPA1 cells. **(A)** Dose-dependent effect of HC030031 on Rb⁺ efflux induced by 50 μM of AITC in the presence of 0.15 and 1.5 mM of Ca²⁺. **(B)** Dose-dependent effect of HC030031 on Rb⁺ efflux induced by 300 μM of AITC in the presence of 0.15, 1.5, and 5 mM of Ca²⁺. **(C)** A plate view of the raw data from a 96-well plate treated with control and 100 μM HC030031. Shown on right are the bar-graphed statistic results for control wells ($n = 42$) and the blocker wells ($n = 42$).

increased AITC to 300 μM and examined the effect of HC030031 in the presence of the different Ca²⁺ concentrations. As shown in *Figure 8B*, in 0.15 mM of Ca²⁺, HC030031 showed a weak effect with the maximal inhibition <10%. Increase in Ca²⁺ to 1.5 mM significantly improved the assay window. The further increase to 5 mM only slightly changed the assay window and the slope of the curve, but elevated both basal and maximal responses. The optimal Ca²⁺ condition for this assay would be 1.5 mM, not only because it provided a better assay window but also because this is within the range of physiological Ca²⁺ concentration. Under this

condition, the IC₅₀ for HC030031 is 30.8 μM. *Figure 8C* shows the raw data from a representative 96-well plate. The plate was treated with the buffer control and 100 μM HC030031 before addition of 300 μM AITC. The percent cellular Rb⁺ contents are 3.0 ± 0.6 and 68.6 ± 2.2, respectively. The Z' factor is calculated to be 0.87, indicating that the assay is robust for HTS.

Taken together, in the present study, we have developed a robust high-throughput Rb⁺ efflux assay for NSCC channels using recombinant TRP channel cell lines. NSCC channels exist in a majority of cell types. The channel permeability to Rb⁺ is sufficient for developing this robust assay. Rb⁺ is often used for detection of potassium transport across membrane. The membrane Rb⁺ fluxes detected by either radioactive or nonradioactive approaches are often interpreted as a reflection of the specific mechanism for K⁺ transportation.^{16,25} Based on the results from this study, NSCC channels are able to carry Rb⁺ upon channel activation. These channels normally provide ionic entries for Na⁺ and Ca²⁺ under physiological conditions. Although they can be equally permeable to K⁺, cells usually have membrane potential close to K⁺ reversal potential when the channels get activated. The K⁺ flux under such conditions is much

weaker than Na⁺ and Ca²⁺. So, the Rb⁺ flux detected using this assay can be a result of the increased NSCC channel activities. Care should be taken when interpreting the data obtained from Rb⁺ efflux assay, especially when working with primary cells.

CONCLUSION

Rb⁺ efflux assay has been demonstrated in this study to be robust for measurement of NSCC channel activities. The examples presented as NSCC channels are TRPC6 and TRPA1, 2 recently identified potential drug targets. Due to the unique activation

properties, the available assays for these channels especially in HTS format are very limited. Rb⁺ efflux assay is sensitive, label-free, and HTS applicable. It provides a robust way to bridge the gap between the extremely low-throughput patch clamp recording and the high-throughput fluorescence-based assay. In addition, it has the potential to be used as an alternative for primary HTS.

ACKNOWLEDGMENTS

We thank Dr. Janet Paulsen and Wei Cao for their help in generating TRPC6 expression construct. We also thank Dr. Ravikumar Peri for providing the HEK293-hTRPA1 stable cell line.

AUTHOR DISCLOSURE STATEMENT

No competing financial interests exist.

REFERENCES

- Abramowitz J, Birnbaumer L: Physiology and pathophysiology of canonical transient receptor potential channels. *FASEB J* 2009;23:297–328.
- Benarroch EE: TRP channels: functions and involvement in neurologic disease. *Neurology* 2008;70:648–652.
- Watanabe H, Murakami M, Ohba T, Ono K, Ito H: The pathological role of transient receptor potential channels in heart disease. *Circ J* 2009;73:419–427.
- Inoue R: TRP channels as a newly emerging non-voltage-gated CA²⁺ entry channel superfamily. *Curr Pharm Des* 2005;11:1899–1914.
- Venkatachalam K, Montell C: TRP channels. *Annu Rev Biochem* 2007;76:387–417.
- Hofmann T, Obukhov AG, Schaefer M, Harteneck C, Gudermann T, Schultz G: Direct activation of human TRPC6 and TRPC3 channels by diacylglycerol. *Nature* 1999;397:259–263.
- Bush EW, Hood DB, Papst PJ, Chapo JA, Minobe W, Bristow MR, et al.: Canonical transient receptor potential channels promote cardiomyocyte hypertrophy through activation of calcineurin signaling. *J Biol Chem* 2006;281:33487–33496.
- Cai R, Ding X, Zhou K, Shi Y, Ge R, Ren G, et al.: Blockade of TRPC6 channels induced G2/M phase arrest and suppressed growth in human gastric cancer cells. *Int J Cancer* 2009;125:2281–2287.
- Guinamard R, Bois P: Involvement of transient receptor potential proteins in cardiac hypertrophy. *Biochim Biophys Acta* 2007;1772:885–894.
- Kuwahara K, Wang Y, McAnally J, Richardson JA, Bassel-Duby R, Hill JA, et al.: TRPC6 fulfills a calcineurin signaling circuit during pathologic cardiac remodeling. *J Clin Invest* 2006;116:3114–3126.
- Nakayama H, Wilkin BJ, Bodi I, Molkentin JD: Calcineurin-dependent cardiomyopathy is activated by TRPC in the adult mouse heart. *FASEB J* 2006;20:1660–1670.
- Onohara N, Nishida M, Inoue R, Kobayashi H, Sumimoto H, Sato Y, et al.: TRPC3 and TRPC6 are essential for angiotensin II-induced cardiac hypertrophy. *EMBO J* 2006;25:5305–5316.
- Shi Y, Ding X, He Z, Zhou K, Wang Q, Wang Y: Critical role of TRPC6 channels in G2 phase transition and the development of human esophageal cancer. *Gut* 2009;58:1443–1450.
- Korsgaard MP, Strøbaek D, Christophersen P: Automated planar electrode electrophysiology in drug discovery: examples of the use of QPatch in basic characterization and high content screening on Na(v), K(Ca)2.3, and K(v)11.1 channels. *Comb Chem High Throughput Screen* 2009;12:51–63.
- Priest BT, Swensen AM, McManus OB: Automated electrophysiology in drug discovery. *Curr Pharm Des* 2007;13:2325–2337.
- Lindström P, Norlund L, Sehlin J: Glucose reduces both Rb⁺ influx and efflux in pancreatic islet cells. *FEBS Lett* 1986;200:67–70.
- Shennan DB, Gow IF: Volume-activated K⁺/Rb⁺ efflux in lactating rat mammary tissue. *Biochim Biophys Acta* 2000;1509:420–428.
- Gill S, Gill R, Wicks D, Liang D: A cell-based Rb⁺-flux assay of the Kv1.3 potassium channel. *Assay Drug Dev Technol* 2007;5:373–380.
- Jow F, Tseng E, Maddox T, Shen R, Kowal D, Dunlop J, et al.: Rb⁺ efflux through functional activation of cardiac KCNQ1/minK channels by the benzodiazepine R-L3 (L-364,373). *Assay Drug Dev Technol* 2006;4:443–450.
- Wang K, McIlvain B, Tseng E, Kowal D, Jow F, Shen R, et al.: Validation of an atomic absorption rubidium ion efflux assay for KCNQ/M-channels using the ion Channel Reader 8000. *Assay Drug Dev Technol* 2004;2:525–534.
- Doerner JF, Gisselmann G, Hatt H, Wetzel CH: Transient receptor potential channel A1 is directly gated by calcium ions. *J Biol Chem* 2007;282:13180–13189.
- Nagata K, Duggan A, Kumar G, Garcia-Añoveros J: Nociceptor and hair cell transducer properties of TRPA1, a channel for pain and hearing. *J Neurosci* 2005;25:4052–4061.
- Zurborg S, Yurgionas B, Jira JA, Caspani O, Heppenstall PA: Direct activation of the ion channel TRPA1 by Ca²⁺. *Nat Neurosci* 2007;10:277–279.
- McNamara CR, Mandel-Brehm J, Bautista DM, Siemens J, Deranian KL, Zhao M, et al.: TRPA1 mediates formalin-induced pain. *Proc Natl Acad Sci USA* 2007;104:13525–13530.
- Gow IF, Thomson J, Davidson J, Shennan DB: The effect of a hyposmotic shock and purinergic agonists on K⁺/Rb⁺ efflux from cultured human breast cancer cells. *Biochim Biophys Acta* 2005;1712:52–61.

Address correspondence to:

Dr. Kun Liu
Department of Screening Sciences
Wyeth Research
500 Arcola Road
Collegeville, PA 19426

E-mail: liuk2@wyeth.com